

Vision Funnels

Product Catalogue

Table of contents




Funnel sizes

Parameters of the different funnel types



Safety funnels with valve

and threaded ports



Safety funnels with valve

for waste canisters

filtered air, pressure release

If the funnel is used without HPLC, the port can be sealed with a blind plug. An exhaust filter must be used to avoid pressure buildup in the container when used in continuous waste disposal mode (e.g.: connected to an HPLC system).



manual valve made of PTFE to provide exeptional chemical resistance

Safety Funnels with Valve*

Funnels with valve for waste canisters with 1/4"-28 open ports

VLF45-0002-TW	DIN45 Safety Cap with safety funnel and 2 ports	1 / package
VLF45-0011-TW	DIN45 Safety Cap with safety funnel and 2 ports	1 / package
VLF61-0002-TW	DIN60/61 Safety Cap with safety funnel and 2 ports	1 / package
VLF61-0011-TW	DIN60/61 Safety Cap with safety funnel and 2 ports	1 / package

* Wetted materials are: PP (Polypropylene), PTFE (Teflon)

harmful vapors

up to 5 capillary

(1/16", 2mm, 3mm,

1/8") for continuous solvent disposal

available in GL45/DIN45 and DIN60/61 size

connections


2x standard

1/4"-28 ports

PTFE body

Safety funnels with valve

for general waste disposal





Manual valve made of PTFE to provide exceptional chemical resistance

Funnels with valve for DIN60/61 canisters		
VLF-X02-61*	DIN60/61 Cap with safety funnel, diameter 155mm with removable sieve	1 / package
VLF-X03-61**	DIN60/61 Cap with safety funnel, diameter 150mm with built-in sieve	1 / package

* Wetted materials are: PP (Polypropylene), PTFE (Teflon), Brass ** Wetted materials are: PP (Polypropylene), PTFE (Teflon)

General laboratory

safety funnels with lid



Safety funnels with lid



Available for different cap sizes: GL45/DIN45, DIN51/S55, DIN60/61, 53B and many more

For available part numbers please check page 14!


PP lid with click-closure



DID YOU KNOW?

With different colours, you can easily identify and store different types of waste such as organic, aqueous, chlorinated, etc.*

* Wetted materials are: PP (Polypropylene), PTFE (Teflon)



burpless adapter inside made of PTFE with excellent chemical resistance consisting SIEVE

harmful vapors are captured inside

Conductive safety funnels



Part No.	Description	Unit
ESD SAFETY F	UNNELS and Canisters	
VLF61-L007	ESD Safety Funnel (black PP) and Lid (black) with DIN60/61 Safety Cap (black)	EA
VLK61-10EL	10L ESD Canister with level sensing (DIN60/61) 300x200x258mm, black, ADR Y1,5	EA
VLK61-20EL	20L ESD Canister with level sensing (DIN60/61) 300x200x495mm, black, ADR Y1,5	EA

VLT-002	Plastic Collecting Tray 300x200x170mm, black
VLT-003	Plastic Collecting Tray 400x300x270mm, black
VLT-004	Plastic Collecting Tray 600x400x170mm, black

ESD COLLECTING TRAYS

Part No.


VLT-001

Description

Plastic Collecting Tray

400x300x170mm, black Plastic Collecting Tray

ESD spill trays







VisionEC-Lab Ltd. | www.visionsafetycap.com

Unit

15I

7,6l

26l 35l

XL funnel with lid

Safety Funnel with Valve*			
Funnel with valve	Funnel with valve for DIN60/61 canisters		
VLF-X01-61	PP Funnel with lid, diameter: 250mm; Including adaptor for DIN60/61 canisters	1 / package	
VLF-XM1-61	PP Funnel with mechanical level sensing and lid, diameter: 250mm; Including adaptor for DIN60/61 canisters	1 / package	

* Wetted materials are: PP (Polypropylene)

CWS - Central Waste Systems


We design and build the waste solvent collection system from easy-to-install elements according to your individual needs.

Integrated waste managing systems are available for laboratory furnitures

Contact us! For more technical information and availability please visit **www.visionsafetycap.com**

OEM requests

are welcome!



Connectors for the CWS system









Three-way tube connector

Union

Elbow connector

Flexible tube connector

CWS adapters





Elbow threaded adapter



Waste container connector

In-line threaded adapter

GL45 threaded adapter



Accessories



Adjustable lenght tube holder



HDPE tubing



Tube holder



Corrugated flexible PP tubing

Notes

..... _____

Part numbers

for Vision Funnels

Safety Funnels v		
	canisters with 1/4"-28 open ports	
VLF45-0002-TW	DIN45 Safety Cap with safety funnel and 2 ports	1 / package
VLF45-0011-TW	DIN45 Safety Cap with safety funnel and 2 ports	1 / package
VLF61-0002-TW	DIN60/61 Safety Cap with safety funnel and 2 ports	1 / package
Funnels with valve for DIN		
VLF-X02-61	DIN60/61 Cap with safety funnel, diameter 155mm with removable sieve	1 / package
VLF-X03-61	DIN60/61 Cap with safety funnel, diameter 150mm with built-in sieve	1 / package
XL Safety Funne		
XL Funnels with Lid for DIM	N60/61 canisters	
VLF-X01-61	PP Funnel with lid, diameter 250mm Including adaptor for DIN60/61 canisters	1 / package
VLF-XM1-61	PP Funnel with mechanical level sensing and lid, diameter 250mm Including adaptor for DIN60/61 canisters	1 / package
Safety Funnels v	vith Lid	
Safety Funnels with Lid for	GL45 bottles	
VLF45-L001	Safety Funnel (neutral PP) and Lid (red) with GL45 Safety Cap (red) Ea	EA
VLF45-L002	Safety Funnel (neutral PP) and Lid (yellow) with GL45 Safety Cap (yellow) Ea	EA
VLF45-L003	Safety Funnel (neutral PP) and Lid (orange) with GL45 Safety Cap (orange) Ea	EA
VLF45-L004	Safety Funnel (neutral PP) and Lid (green) with GL45 Safety Cap (green) Ea	EA
VLF45-L005	Safety Funnel (black PP) and Lid (black) with GL45 Safety Cap (black) Ea	EA
VLF45-L006	Safety Funnel (black PP) and Lid (black) with GL45 Safety Cap (yellow) Ea	EA
Safety Funnels with Lid for	DIN45 canisters	
VLF45-L007	Safety Funnel (neutral PP) and Lid (red) with DIN45 Safety Cap (red)	EA
VLF45-L008	Safety Funnel (neutral PP) and Lid (yellow) with DIN45 Safety Cap (yellow)	EA
VLF45-L009	Safety Funnel (neutral PP) and Lid (orange) with DIN45 Safety Cap (orange)	EA
VLF45-L010	Safety Funnel (neutral PP) and Lid (green) with DIN45 Safety Cap (green)	EA
VLF45-L011	Safety Funnel (black PP) and Lid (black) with DIN45 Safety Cap (black)	EA
VLF45-L012	Safety Funnel (black PP) and Lid (black) with DIN45 Safety Cap (yellow)	EA
Safety Funnels with Lid for	53B canisters	
VLF53-L001	Safety Funnel (neutral PP) and Lid (red) with 53B Safety Cap (red)	EA
VLF53-L002	Safety Funnel (neutral PP) and Lid (yellow) with 53B Safety Cap (yellow)	EA
VLF53-L003	Safety Funnel (neutral PP) and Lid (orange) with 53B Safety Cap (orange)	EA
VLF53-L004	Safety Funnel (neutral PP) and Lid (green) with 53B Safety Cap (green)	EA
VLF53-L005	Safety Funnel (black PP) and Lid (black) with 53B Safety Cap (black)	EA
VLF53-L006	Safety Funnel (black PP) and Lid (black) with 53B Safety Cap (yellow)	EA
Safety Funnels with Lid for	DIN51/S55 canisters	
VLF55-L001	Safety Funnel (neutral PP) and Lid (red) with DIN51/S55 Safety Cap (red)	EA
VLF55-L002	Safety Funnel (neutral PP) and Lid (yellow) with DIN51/S55 Safety Cap (yellow)	EA
VLF55-L003	Safety Funnel (neutral PP) and Lid (orange) with DIN51/S55 Safety Cap (orange)	EA
VLF55-L004	Safety Funnel (neutral PP) and Lid (green) with DIN51/S55 Safety Cap (green)	EA
VLF55-L005	Safety Funnel (black PP) and Lid (black) with DIN51/S55 Safety Cap (black)	EA
Safety Funnels with Lid for	DIN60/61 canisters	
VLF61-L001	Safety Funnel (neutral PP) and Lid (red) with DIN60/61 Safety Cap (red)	EA
VLF61-L002	Safety Funnel (neutral PP) and Lid (yellow) with DIN60/61 Safety Cap (yellow)	EA
VLF61-L003	Safety Funnel (neutral PP) and Lid (orange) with DIN60/61 Safety Cap (orange)	EA
VLF61-L004	Safety Funnel (neutral PP) and Lid (green) with DIN60/61 Safety Cap (green)	EA
VLF61-L005	Safety Funnel (black PP) and Lid (black) with DIN60/61 Safety Cap (black)	EA
VLF61-L006	Safety Funnel (black PP) and Lid (black) with DIN60/61 Safety Cap (yellow)	EA
VLF61-L007	ESD Safety Funnel (black PP) and Lid (black) with DIN60/61 Safety Cap (black)	EA

For more information visit www.visionsafetycap.com

Contact email: info@visionec.hu

HUNGARY, H-2132 Göd, Attila str. 13.

